

GEL MOBILITY SHIFT ASSAY

Material/reagents:

- 2X binding buffer:
20 mM Tris pH 7.5
100 mM NaCl
2 mM EDTA
10% glycerol
- poly(dI-dC) (5 μ g/ μ l in TE)
- BSA (10 mg/ml)
- 32 P-labeled probe
- Buffer C:
20 mM HEPES pH 7.9
25% glycerol
420 mM NaCl
1.5 mM $MgCl_2$
0.5 mM DTT
0.5 mM PMSF
0.2 mM EDTA
- 4% native polyacrylamide gel (50 ml):
5 ml 40% acrylamide (30:1)
2.5 ml 5X TBE
41.95 ml H_2O
0.5 ml 10% APS
50 μ l TEMED
- 0.25X TBE running buffer
- Loading dye:
0.2% bromophenol blue
0.2% xylene cyanol
50% glycerol

Procedure:

1. Prepare reaction mix (20 μ l/sample):

dH ₂ O	6.3 μ l
2x binding buffer	12.5 μ l
BSA (10 mg/ml)	1.0 μ l
poly(dI-dC) (5 μ g/ μ l)	0.1 μ l
probe (1 ng/ μ l)	0.1 μ l

2. Add buffer C to the reaction mix if your volume of whole cell extract will be less than 5 μ l. The sum of whole cell extract and buffer C should be 5 μ l.

3. Add 10–15 μ g of whole cell extract ($\leq$ 5 μ l) to the reaction mix and resuspend. The total volume should be 25 μ l.

4. Incubate at room temperature for 20 min.

5. Add 2.5 μ l of loading dye.

6. Load sample onto a 4% polyacrylamide gel.

7. Run the gel at room temperature for ~2.5 hours at 120-130 V in 0.25X TBE. (Run until the bromophenol blue is ~1 inch from the bottom of the gel. Free probe will co-migrate with the blue dye.)

8. Dry the gel at 80°C and expose to film or PhosphoImager.

Comments:

This protocol works well for HSF DNA-binding activity experiments. An alternative protocol by Gen Matsumoto uses a 19:1 40% polyacrylamide and no BSA, but is otherwise the same.

References:

1. Mosser et al., 1988, Mol. Cell. Biol. 8:4736-4744.
2. Fried and Crothers, 1981, Nucleic Acids Res. 9:6505-6525
3. Garner and Revzin, 1981, Nucleic Acids Res. 9:3047-3060
4. Singh et al., 1986, Nature 319:154-158

Submitted by:

Carina Holmberg